

SPE Cartridges



InnoSep

New

Please contact Life Science
for more information
888-226-2775

- **Easy to use** - sample preparation in four simple steps
- **Fast** - nine times faster than traditional extraction techniques
- **Economical, environmentally friendly**
 - reduces use of hazardous solvents by 95%
 - reduces operator time, glassware and equipment by 83%
- **Reliable**
 - InnoSep offers the tightest quality control for particle size, surface activity, pH, flow rate and interferences.
 - hermetically sealed packets prevent contamination and allows inspection of the cartridges before use.

Ordering Information



¹ Encapped

Quantity	100/box	50/box	50/box	30/box	30/box
Phase	100mg 1mL	200mg 3mL	500mg 3mL	500mg 6mL	1000mg 6mL
C18 Octadecyl C18/14%	IS72001	IS72002	IS72003	IS72006	IS72007
C18 Octadecyl C18/18%	IS12001	IS12002	IS12003	IS12006	IS12007
C18 Octadecyl C18/22%	IS12101	IS12102	IS12103	IS12106	IS12107
C18 Octadecyl C18/OH	IS12201	IS12202	IS12203	IS12206	IS12207
C8 Octyl	IS72011	IS72012	IS72013	IS72016	IS72017
PH Phenyl	IS72081	IS72082	IS72083	IS72086	IS72087
CH Cyclohexyl	IS72071	IS72072	IS72073	IS72076	IS72077
C4 Butyl	IS72021	IS72022	IS72023	IS72026	IS72027
C2 Ethyl	IS72061	IS72062	IS72063	IS72066	IS72067
C1 Methyl	IS72051	IS72052	IS72053	IS72056	IS72057
Price	\$155.00	\$116.00	\$130.00	\$ 86.00	\$122.00
SI Silica Gel	IS72101	IS72102	IS72103	IS72106	IS72107
FLO Florisil	IS72111	IS72112	IS72113	IS72116	IS72117
ALN Alumina (Neutral)	IS72121	IS72122	IS72123	IS72126	IS72127
ALA Alumina (Acidic)	IS72131	IS72132	IS72133	IS72136	IS72137
ALB Alumina (Basic)	IS72141	IS72142	IS72143	IS72146	IS72147
Price	\$150.35	\$ 99.20	\$116.25	\$ 79.00	\$102.30
CNe Cyanopropyl (ec) ¹	IS72221	IS72222	IS72223	IS72226	IS72227
CN Cyanopropyl	IS72201	IS72202	IS72203	IS72206	IS72207
DIO Diol	IS72151	IS72152	IS72153	IS72156	IS72157
COOH Carboxylic Acid	IS72311	IS72312	IS72313	IS72316	IS72317
SCX Benzenesulfonic Acid	IS72321	IS72322	IS72323	IS72326	IS72327
NH ₂ Aminopropyl	IS72211	IS72212	IS72213	IS72216	IS72217
PSA Pri/Sec Amino	IS72241	IS72242	IS72243	IS72246	IS72247
DEA Diethylamino	IS72331	IS72332	IS72333	IS72336	IS72337
N+ Quaternary Amino	IS72301	IS72302	IS72303	IS72306	IS72307
Price	\$212.35	\$139.50	\$170.50	\$114.70	\$158.10

SPE 96-Well Plates

InnoSep 96-Well Plates



Solid Phase Extraction in a 96-Well Format

•Economic

- prepare 96 samples at one time
- reduce use of hazardous solvents, operator time, glassware and equipment.

•Reproducible


- achieve high recoveries with 2-3% coefficients of variation

•Versatile

- perform extractions with a wide variety of sorbents.
- two plate sizes and a range of standard packing amounts are available.

Ordering Information

Phase	25mg 1mL	50mg 1mL	25mg 2mL	50mg 2mL	100mg 2mL	200mg 2mL
C18 Octadecyl C18/14%	IS24050	IS24051	IS24052	IS24053	IS24054	IS24055
C18 Octadecyl C18/18%	IS24060	IS24061	IS24062	IS24063	IS24064	IS24065
C18 Octadecyl C18/22%	IS24070	IS24071	IS24072	IS24073	IS24074	IS24075
C18 Octadecyl C18/OH	IS24080	IS24081	IS24082	IS24083	IS24084	IS24085
C8 Octyl	IS24100	IS24101	IS24102	IS24103	IS24104	IS24105
PH Phenyl	IS24110	IS24111	IS24112	IS24113	IS24114	IS24115
CH Cyclohexyl	IS24120	IS24121	IS24122	IS24123	IS24124	IS24125
C4 Butyl	IS24130	IS24131	IS24132	IS24133	IS24134	IS24135
C2 Ethyl	IS24140	IS24141	IS24142	IS24143	IS24144	IS24145
C1 Methyl	IS24150	IS24151	IS24152	IS24153	IS24154	IS24155
CNe Cyanopropyl (ec) ¹	IS24160	IS24161	IS24162	IS24163	IS24164	IS24165
CN Cyanopropyl	IS24170	IS24171	IS24172	IS24173	IS24174	IS24175
SI Silica Gel	IS24180	IS24181	IS24182	IS24183	IS24184	IS24185
FLO Florisil	IS24190	IS24191	IS24192	IS24193	IS24194	IS24195
ALN Alumina (Neutral)	IS24200	IS24201	IS24202	IS24203	IS24204	IS24205
ALA Alumina (Acidic)	IS24210	IS24211	IS24212	IS24213	IS24214	IS24215
ALB Alumina (Basic)	IS24220	IS24221	IS24222	IS24223	IS24224	IS24225
DIO Diol	IS24230	IS24231	IS24232	IS24233	IS24234	IS24235
COOH Carboxylic Acid	IS24240	IS24241	IS24242	IS24243	IS24244	IS24245
SCX Benzenesulfonic Acid	IS24250	IS24251	IS24252	IS24253	IS24254	IS24255
NH ₂ Aminopropyl	IS24260	IS24261	IS24262	IS24263	IS24264	IS24265
PSA Pri/Sec Amino	IS24270	IS24271	IS24272	IS24273	IS24274	IS24275
DEA Diethylamino	IS24280	IS24281	IS24282	IS24283	IS24284	IS24285
N+ Quaternary Amino	IS24290	IS24291	IS24292	IS24293	IS24294	IS24295
Price	\$ 286.75	\$ 317.75	\$ 302.25	\$ 333.25	\$ 348.75	\$ 364.25



1. Condition by aspirating a few microliters of suitable solvent or solution through the column.
2. Extract and concentrate compounds of interest by aspirating a sample solution.
3. Remove co-extracted interferences by aspirating a few microliters of suitable solvent or solution.
4. Elute extracted compounds by aspirating a few microliters of suitable solvent or solution.

¹ Endcapped